



Send an email to “[thomas.deegan22@imperial.ac.uk](mailto:thomas.deegan22@imperial.ac.uk)”  
with your github username.

Thomas will add you to the competition:

- Google Drive shared folder
- Github repository



David Wang  
[david.wang@elemindtech.com](mailto:david.wang@elemindtech.com)

June 26th, 2025



### Action Button Tips

Double-Tap starts a pre-programmed routine, LED will glow purple for 5 seconds.  
Send command "ble\_therapy 0" to stop.

### LED Light Tips

Blue Breathing = BLE Pairing To Phone Needed  
Red+Blue Alternating = Needs Real-World Time  
Green Breathing = Charging

Press any button when off charger to check charge:  
Green = >70%+ Charged  
Orange = 25-70% Charged  
Red = <25% Charged

# How to Wear



## Tips:


- Wipe forehead and over ears with napkin to remove oils.
- Does not need to be tight. Just not shifting if you shake your head.
- For long hair:
  - Set adjustable strap to largest size.
  - Slip over hair down to neck.
  - Bring headband up to the forehead, under hair.
  - Tighten strap as need.
- Once on the head:
  - Clear hair from between references and ear.
  - Pull down on the adjustable strap to ensure references are well-seated on the ears.

## Hardware Specs


- **Dry EEG Sensors:** 3-channels (FP1, FP2, FPZ), 2 reference electrodes  
(Bandwidth: 0.5-35Hz, Sample Rate: 250 Hz)
- **3-axis Accelerometer** (+/-2Gs, Sample Rate: 25 Hz)
- **Stimulation:** Phase-locked acoustic stimulation via bone conduction transducer
- **Real-time phase computation, stimulation, and ML sleep-staging (on-device)**
- **Phase Stimulation Accuracy:** 90-98% PLV (from peer-reviewed study)
- **Phase Stimulation Error:** 4-15 degrees (from peer-reviewed study)
- **EEG-based sleep tracking accuracy:** 83%
- **Weight:** 60g
- **Battery Life:** >3hr w/ stim, >24 hrs recording, 48hr standby time
- **Battery Charge Time:** ~3 hrs
- **Connectivity:** Bluetooth Connectivity, USB-C
- **Automatic Data Transfer:** Lossy or Lossless sensor data transferred over BLE to App to Cloud.



# Closed Loop Audio Stimulation (CLAS)



# Command Overview




# Morpheus Controller (A Java GUI for Testing)

File  
Upload/  
Download  
Menu

Connect to Headband

Local Logging  
(Status text saved to  
Downloads folder  
automatically)


One-Click  
Convenience  
Commands





# Morpheus Controller

## (A Java GUI for Testing)



1. Select the port  
Look for “Elemind Morpheus”.  
On Mac use either cu.usbmodem or tty.usbmodem.
2. Click Open
3. Type Commands (1 command per line)
4. Click “Send Cmds”

# Morpheus Controller - Impedance Check

The screenshot shows the Morpheus Controller v1.80.ext application window. The interface includes a menu bar (Configuration, Window, Help), a Port selection field (tty.usbmodem1101 (Elemind Morpheus (Dial-In))), and a Log To field. The main display area has tabs for EEG, LeadOff, Stimulation, Accelerometer, Audio, Sleep Stage, and Status. The LeadOff tab is active, showing two plots: 'LeadOff' (raw value vs Time [s]) and 'LeadOff Inter Sample Period' (count vs time [us]). The 'LeadOff' plot shows a green line (fp2) at approximately 650, a purple line (fp1) at approximately 600, and a yellow line (ac\_excite) at approximately 100. The 'LeadOff Inter Sample Period' plot shows a blue line (mode) at approximately 10. The 'Commands' panel on the right contains buttons for Test EEG, Skip BLE Bond, Set Time, Load Cmds, Save Cmds, and Send Cmds. Arrows point from the numbered instructions to the corresponding UI elements.

1. Select the port  
Look for “Elemind Morpheus”.  
On Mac use either cu.usbmodem or tty.usbmodem.
2. Click Open
3. Click “Test EEG”
4. Click LeadOff Tab

Leadoff values range from 0~22000  
A value below 16000 is okay.

# Morpheus Controller Setup

## First Time Setup


- Install Java Development Kit (JDK) 24 from Oracle

(<https://www.oracle.com/uk/java/technologies/downloads/>)

It has better USB compatibility and streams data faster.


- Mac Security:

- The first time you run Morpheus Controller, you will see this. See next slide for fix.



# Morpheus Controller Setup

- Mac Security (Continued):
  - To Allow, Navigate to: Apple Icon (Upper Left) -> System Settings -> Security & Privacy



- Scroll Down and click “Open Anyway”

